


Digital Plethysmometer

Model No. VJDP-01



Foot Switch

TECHNICAL SPECIFICATIONS



Video Link

General

The Digital Plethysmometer is useful Apparatus to measure the changes in small paw volume indicating inflammation. The Apparatus measures micro volume changes due to displacement of water with the insertion of Rat or Mouse paw.

Operation

- In Comparison with other digital Plethysmometer, we use advance imaging technology with Artificial Intelligence (AI).
- Simple User Interface
- Water loss correction with each reading
- Image & AI Based technology.
- Foot Switch (USB)
- Software
- Accuracy - 0.01 ml
- 35mm Tube diameter for Mouse & Rat
- No need of Sensors (it will avoid cost of replacement of sensor)

Technical Specifications.

Sr. No.	Specification	Details
1	Certificate	GMP, ISO, CE
2	Material of Construction	SS 316, Glass , Acrylic
3	Tube diameter	35 mm diameter (Glass)
4	Calibration Cuber	1 CC
5	System Requirement	OS Windows 8 and above
6	Resolution	0.01 ml
7	Dimension (LXWXH) in mm	425mm X 280 mm X 300 mm approximately
8	Weight	11 kg 35 kg with packing.

Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary