

# Mazes

**Light / Dark Box  
Place Preference Box**

## TECHNICAL SPECIFICATIONS

## Light Dark Box




Sr. No.	Details
1	Available For Rat & Mouse
2	Dimension - Mouse, 20cm x 40cm (each side). Rat, 50cm x 50cm (each side)
3	Base- Grey / Black (Non - Reflective)
4	Material - Acrylic
5	Model No. -For Mouse- VJLDB-01M For Rat- VJLDB-01R

### **Note -**

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary

## Place Preference Box



Sr. No.	Details
1	Available For Rat & Mouse
2	Overall dimensions: Mouse 40cm x 40cm.
3	Overall dimensions: Rat 100cm x 50cm
4	Available with transparent or opaque walls
5	Grey / Black (Non - Reflective)
6	Material - Acrylic / Aluminium
7	Model No. for Mouse -VJPPB-01M / for Rat-VJPPB-01R

### Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary