

Mazes

Radial Maze
Elevated Plus Maze
Zero Maze
Y-Maze
T-Maze

TECHNICAL SPECIFICATIONS

Radial Maze



Sr. No.	Details
1	Available For Rat & Mouse
2	Rat Maze ARM Dimension 50 cm X10 cm (Wall Height 20 cm)
3	Mouse Maze ARM Dimension 35 cm X 5 cm (Wall Height 10 cm)
4	8 Arm
5	Grey / Black (Non - Reflective)
6	Material - Mild Steel / Aluminium
7	Model No.For Mouse - VJRM-01M For Rat- VJRM-01R

Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary

Elevated Plus Maze



Sr. No.	Details
1	Available For Rat & Mouse
2	Elevated 50 cm from Floor
3	Rat Maze ARM Dimension 50 cm X10 cm (Wall Height Close Arm 20 cm)
4	Mouse Maze ARM Dimension 35 cm X 5 cm (Wall Height Close Arm 10 cm)
5	Grey / Black (Non - Reflective)
6	Material - Mild Steel / Aluminium
7	Model No. for Mouse -VJEPM-01M / for Rat-VJEPM-01R

Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary

Zero Maze



Sr. No.	Details
1	Available For Rat & Mouse
2	Elevated 50 cm from floor
3	Rat maze dimension 100 cm diameter, 10 cm lane width, 20 cm wall height
4	Mice maze dimension 50 cm diameter, 5 cm lane width, 15 cm wall height
5	Infrared transparent inner walls which allow a single camera to see the animal in all parts of the maze.
6	Grey / Black (Non - Reflective)
7	Material - Mild Steel / Aluminium
8	Model No. for Mouse -VJZM-01M/ for Rat-VJZM-01R

Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary

Y-Maze



Sr. No.	Details
1	Available For Rat & Mouse
2	Rat Maze ARM Dimension 50 cm X10 cm (Wall Height 20 cm)
3	Mouse Maze ARM Dimension 35 cm X 5 cm (Wall Height 10 cm)
4	3 Arm
5	Grey / Black (Non - Reflective)
6	Material - Mild Steel / Aluminium
7	Model No. for Mouse -VJYM-01M/ for Rat-VJYM-01R

Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary

T-Maze



Sr. No.	Details
1	Available For Rat & Mouse
2	Rat Maze ARM Dimension 50 cm X10 cm (Wall Height 20 cm)
3	Mouse Maze ARM Dimension 35 cm X 5 cm (Wall Height 10 cm)
4	3 Arm
5	Grey / Black (Non - Reflective)
6	Material - Mild Steel / Aluminium
7	Model No. for Mouse -VJTM-01M/ for Rat-VJTM-01R

Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary