

# Mazes

**Morris Water Maze**  
**Barnes Maze**

## TECHNICAL SPECIFICATIONS

## Morris Water Maze



Sr. No.	Details
1	Available For Rat & Mouse
2	Dimension - Rat 180cm, Rat/Mouse 150cm, Mouse 120cm
3	Casters with wheel brakes
4	Elevated with drainage system
5	Islands with Hight adjustment for use in Morris Maze applications
6	Material - FRP (Fibre rain-forced plastic)
7	Model No.For Mouse- VJMWM-01M For Rat- VJMWM-01R For Mouse/Rat - VJMWM-MR-150

### Note -

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary

## **Barnes Maze**

---

Sr. No.	Details
1	Available For Rat & Mouse
2	Elevated to convenient working height.
3	Dimensions: Rat 122cm, Mouse 91cm
4	Holes - 20
5	Grey / Black (Non - Reflective)
6	Material - Mild Steel / Aluminium
7	Model No. for Mice -VJEPM-01M / for Rat-VJEPM-01R

### **Note -**

1. All the images displayed in the offer are for representation purpose only but actual may vary.
2. Dimension & weight may vary